d Envi-...itatia **CA56 Cylinder Analysis System** . planketing ustion Analysing

ciency Ultraviolet

ty Gas Production

Powdered Metals lled Environments

ens Manufacturing

rocarbon Refining ■ Food Packaging cics ■ Glove Boxes

tron Beam ■ R & D s Fermentation

led Environments ns Manufacturing

Vessel Blanketing

ustion Analysing

iency Ultraviolet y Gas Production

owdered Metals d Environments

Manufacturing

isel Blanketing

ion Analysing

cy Ultraviolet as Production


dered Metals nvironments

anufacturing

l Blanketing



Cylinder Analysis System for measurement of trace gas impurities in inert gas cylinders



Features & Benefits

- Simple high pressure (250 Bar) connection to the cylinder
- Quickly determine trace gas impurities within minutes.
- Can be used with any gas analyzer

Technical Specifications

Test Cylinder Pressure

Gas Fittings

Enclosure

Dimensions

Weight

2-250 barg (15-3700 PSI) regulated

1/8 in Swagelok® (supplied)

Epoxy coated heavy gauge steel

9.3W x 15.4D x 7.5H in (237W x 390D x 190H mm)

15.4 lb (7 kg)

The cylinder analysis system from Systech Illinois is capable of quickly and easily measuring trace gas impurities such as oxygen, moisture or other trace impurities. This is essential for inert gas manufacturers who require every gas cylinder to be checked for impurities before delivery to their customers.

The issue with using trace oxygen and moisture analysers in this type of application is the time

required for the cells to purge down to low levels. To purge the system from ambient air can take several hours or even days. This makes the analysis of many cylinders per day very challenging.

Therefore Systech Illinois have developed a cylinder purge sample system to make these connections simple, to purge each cylinder and to determine the trace impurities within minutes.